

Help Button with Call Acknowledge

Features & Benefits

- Press anywhere on the face for easy alarm activation
- Large Acknowledge button makes it easy to record when staff responds to a call
- Easy to clean membrane switch
- Supervised

Specifications

Part Number:	1840249
Color:	White
Face:	Grey "Alert" face with discreet "Acknowledge" button across top
Dimensions:	2.56" wide, 4" high, .94" deep
Batteries:	3.6 VDC Lithium
Average Battery Life:	5-8 years
Compatibility:	All Lifeline Wireless E-Call Learn Mode Systems



Lifeline Systems
111 Lawrence Street
Framingham, MA 01702
Telephone: 800-451-0525
Fax: 508-988-1384
www.lifelinesystems.com/seniorliving

Lifeline Systems Canada, Inc.
95 Barber Greene Road, Suite 105
Toronto, ON M3C 3E9, Canada
Telephone: 800-387-8120
Fax: 416-445-5402
www.lifelinecanada.com