CAREPOINT

Passive Infared (PIR) Motion Sensor

Assure resident safety and efficient staff operation with passive resident check-in. Central Monitoring Station checks for resident activity during specified period and notifies staff if resident is in-active.

Features & Benefits

- Can be utilized for automated resident check-in or activity monitoring
- A wall tamper triggers an alarm if the sensor is removed
- 3-minute lockout timer extends battery life
- Wall-mounting bracket
- Supervised

Specifications

Color: White

Dimensions: 2.75" wide, 4.38" high, 1.75" deep

Battery: 3.5 VDC lithium battery

Average battery life: 4-5 years
Installation: Wall mounting

Corporate Offices

111 Lawrence Street Framingham,MA 01702 Telephone: 800-816-4885 Fax: 508-988-1363 www.lifelineseniorliving.com 95 Barber Greene Road, Suite 105 Toronto, ON M3C 3E9, Canada Telephone: 800-387-8120 Fax: 416-445-5402 www.lifeline.ca

Service and Support Offices

111 Lawrence Street Framingham,MA 01702 3975-A Industrial Way Concord, CA 94520

1325 Northmeadow Parkway Roswell, GA 30076



MKT 1830, Rev. 03

